

BS EN 62464-2:2011



BSI Standards Publication

Magnetic resonance equipment for medical imaging

Part 2: Classification criteria for
pulse sequences

bsi.

...making excellence a habit.™

National foreword

This British Standard is the UK implementation of EN 62464-2:2011. It is identical to IEC 62464-2:2010.

The UK participation in its preparation was entrusted by Technical Committee CH/62, Electrical Equipment in Medical Practice, to Subcommittee CH/62/2, Diagnostic imaging equipment.

A list of organizations represented on this committee can be obtained on request to its secretary.

This publication does not purport to include all the necessary provisions of a contract. Users are responsible for its correct application.

© BSI 2011

ISBN 978 0 580 69735 7

ICS 11.040.55

Compliance with a British Standard cannot confer immunity from legal obligations.

This British Standard was published under the authority of the Standards Policy and Strategy Committee on 31 March 2011.

Amendments issued since publication

Amd. No.	Date	Text affected
-----------------	-------------	----------------------

**Magnetic resonance equipment for medical imaging -
Part 2: Classification criteria for pulse sequences
(IEC 62464-2:2010)**

Appareils à résonance magnétique utilisés
pour l'imagerie médicale -
Partie 2: Critères de classification pour les
séquences d'impulsions
(CEI 62464-2:2010)

Magnetresonanzgeräte für die
medizinische Bildgebung -
Teil 2: Klassifizierungskriterien für
Pulssequenzen
(IEC 62464-2:2010)

This European Standard was approved by CENELEC on 2011-02-01. CENELEC members are bound to comply with the CEN/CENELEC Internal Regulations which stipulate the conditions for giving this European Standard the status of a national standard without any alteration.

Up-to-date lists and bibliographical references concerning such national standards may be obtained on application to the Central Secretariat or to any CENELEC member.

This European Standard exists in three official versions (English, French, German). A version in any other language made by translation under the responsibility of a CENELEC member into its own language and notified to the Central Secretariat has the same status as the official versions.

CENELEC members are the national electrotechnical committees of Austria, Belgium, Bulgaria, Croatia, Cyprus, the Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, the Netherlands, Norway, Poland, Portugal, Romania, Slovakia, Slovenia, Spain, Sweden, Switzerland and the United Kingdom.

CENELEC

European Committee for Electrotechnical Standardization
Comité Européen de Normalisation Electrotechnique
Europäisches Komitee für Elektrotechnische Normung

Management Centre: Avenue Marnix 17, B - 1000 Brussels

Foreword

The text of document 62B/807/FDIS, future edition 1 of IEC 62464-2, prepared by SC 62B, Diagnostic imaging equipment, of IEC TC 62, Electrical equipment in medical practice, was submitted to the IEC-CENELEC parallel vote and was approved by CENELEC as EN 62464-2 on 2011-02-01.

Attention is drawn to the possibility that some of the elements of this document may be the subject of patent rights. CEN and CENELEC shall not be held responsible for identifying any or all such patent rights.

The following dates were fixed:

- latest date by which the EN has to be implemented at national level by publication of an identical national standard or by endorsement (dop) 2011-11-01
- latest date by which the national standards conflicting with the EN have to be withdrawn (dow) 2014-02-01

In this standard, the following print types are used:

- Requirements and definitions: roman type.
- *Test specifications: italic type.*
- Informative material appearing outside of tables, such as notes, examples and references: in smaller type. Normative text of tables is also in a smaller type.
- TERMS DEFINED IN CLAUSE 3 OF THE GENERAL STANDARD, IN THIS PARTICULAR STANDARD OR AS NOTED: SMALL CAPITALS.

The verbal forms used in this standard conform to usage described in Annex H of the ISO/IEC Directives, Part 2. For the purposes of this standard, the auxiliary verb:

- “shall” means that compliance with a requirement or a test is mandatory for compliance with this standard;
- “should” means that compliance with a requirement or a test is recommended but is not mandatory for compliance with this standard;
- “may” is used to describe a permissible way to achieve compliance with a requirement or test.

Annex ZA has been added by CENELEC.

Endorsement notice

The text of the International Standard IEC 62464-2:2010 was approved by CENELEC as a European Standard without any modification.

Annex ZA (normative)

Normative references to international publications with their corresponding European publications

The following referenced documents are indispensable for the application of this document. For dated references, only the edition cited applies. For undated references, the latest edition of the referenced document (including any amendments) applies.

NOTE When an international publication has been modified by common modifications, indicated by (mod), the relevant EN/HD applies.

<u>Publication</u>	<u>Year</u>	<u>Title</u>	<u>EN/HD</u>	<u>Year</u>
IEC 60601-2-33	2010	Medical electrical equipment - Part 2-33: Particular requirements for the basic safety and essential performance of magnetic resonance equipment for medical diagnosis	EN 60601-2-33 + corr. October	2010 2010
IEC/TR 60788	2004	Medical electrical equipment - Glossary of defined terms	-	-

CONTENTS

INTRODUCTION 5

1 Scope..... 6

2 Normative references..... 6

3 Terms and definitions..... 6

4 PULSE SEQUENCE classification 7

 4.1 General..... 7

 4.2 PULSE SEQUENCE type 7

 4.2.1 General..... 7

 4.2.2 Notation 8

 4.3 Magnetisation modification..... 8

 4.3.1 General..... 8

 4.3.2 Notation 8

 4.4 Dimensionality 10

 4.4.1 General..... 10

 4.4.2 Notation 10

 4.5 Echo number 10

 4.5.1 General..... 10

 4.5.2 Notation 10

Annex A (informative) Examples of use of the PULSE SEQUENCE classification..... 11

Bibliography..... 12

Index of defined terms used in this standard 13

Table 1 – Magnetisation modification techniques 9

Table A.1 – MANUFACTURER-specific classification examples 11

INTRODUCTION

Presently the MANUFACTURERS of MR EQUIPMENT use names for PULSE SEQUENCES which are adopted from the literature (e.g. SPIN-ECHO) or are defined by the MANUFACTURER (e.g. FISP: fast imaging with steady state precession). In the absence of a classification standard for PULSE SEQUENCES, the MANUFACTURER-specific terminology complicates comparison of PULSE SEQUENCES.

The DICOM standard allows the inclusion of PULSE SEQUENCE information with digital MAGNETIC RESONANCE (MR) images. This information helps with the interpretation of images. However, the DICOM standard allows MANUFACTURER-specific terminology.

This International Standard specifies a concise MANUFACTURER-independent classification scheme for MR imaging PULSE SEQUENCES.

In terms of MR imaging, the PULSE SEQUENCE is a chronology of RF-pulses, switching of gradient fields and data acquisition with the intention to create one or more images. As the exact chronology determines the image contrast, image artefacts and other properties of the image, it is necessary to define a consistent and accurate PULSE SEQUENCE classification.

The proposed PULSE SEQUENCE classification notation could be implemented as a new DICOM tag in addition to the existing MANUFACTURER-specific PULSE SEQUENCE name. This would facilitate end users' access to this information. Implementation as a new tag would ensure backward compatibility.

MAGNETIC RESONANCE EQUIPMENT FOR MEDICAL IMAGING –

Part 2 – Classification criteria for pulse sequences

1 Scope

This International Standard specifies the description of PULSE SEQUENCES of MAGNETIC RESONANCE imaging.

NOTE The classification in this standard is suitable for:

- tender texts;
- image annotation;
- protocol definition;
- technical publications.

This International Standard does not apply to MAGNETIC RESONANCE spectroscopy. The classification does not focus on image contrast (T1, T2, proton density), as this is defined by PULSE SEQUENCE parameters (e.g. repetition time, echo time) and is not a property of the PULSE SEQUENCE alone. The PULSE SEQUENCE classification does not specify the K-SPACE acquisition scheme, reconstruction algorithm or post-processing.

2 Normative references

The following referenced documents are indispensable for the application of this document. For dated references, only the edition cited applies. For undated references, the latest edition of the referenced document (including any amendments) applies.

IEC 60601-2-33:2010, *Medical electrical equipment – Part 2-33: Particular requirements for the basic safety and essential performance of magnetic resonance equipment for medical diagnosis*

IEC 60788:2004, *Medical electrical equipment – Glossary of defined terms*

3 Terms and definitions

For the purposes of this document, the terms and definitions given in IEC 60601-2-33:2010, IEC 60788:2004 and the following apply.

3.1

pulse sequence

chronology of radiofrequency-pulses, switching of magnetic field gradients, and data acquisition for the generation of one or more MAGNETIC RESONANCE images

NOTE The terms “imaging sequence” or “sequence” are sometimes used as synonyms for PULSE SEQUENCE.

3.2

transverse magnetisation

magnetisation component perpendicular to the direction of the static magnetic field

3.3

longitudinal magnetisation

magnetisation component parallel to the direction of the static magnetic field

3.4

k-space

mathematical space in which the Fourier transform of the image data is represented

NOTE This is partially or totally filled with the sampled measurement data.

3.5

spin-echo

SE

refocused TRANSVERSE MAGNETISATION arising at time T after an excitation RF pulse and an additional RF-pulse at time $T/2$

NOTE The excitation pulse is typically a 90° RF pulse, and the additional refocusing pulse is typically a 180° RF pulse. SE's can be refocused using a sequence of additional RF pulses.

3.6

gradient-echo

GR

refocused TRANSVERSE MAGNETISATION after a RF pulse using re-phasing magnetic field gradients

4 PULSE SEQUENCE classification

4.1 General

In an imaging PULSE SEQUENCE the LONGITUDINAL MAGNETISATION is partially or totally converted into TRANSVERSE MAGNETISATION via a radio-frequency (RF) excitation pulse. The TRANSVERSE MAGNETISATION is precessing with the Larmor frequency. The precessing TRANSVERSE MAGNETISATION induces the MR signal. For spatial encoding the precessing TRANSVERSE MAGNETISATION is phase-encoded with the help of gradients prior to the data acquisition, and frequency encoding is utilised during the data acquisition. The acquired signal is then stored in a so-called K-SPACE line of the raw data matrix.

The K-SPACE can be two- (2D) or three-(3D)-dimensional. There are several algorithms that allow reconstructing images from incomplete K-SPACE data sets (half or partial Fourier, parallel imaging techniques) – for the PULSE SEQUENCE classification these techniques are not considered.

PULSE SEQUENCES are categorised according to the following classifiers:

- a) magnetisation modification (optional): modification(s) of the LONGITUDINAL MAGNETISATION OR TRANSVERSE MAGNETISATION;
- b) PULSE SEQUENCE type: SPIN-ECHO or GRADIENT-ECHO with the number repetitions (or shots) and the number of K-SPACE lines per RF excitation;
- c) dimensionality of data acquisition: 2D or 3D;
- d) echo number (optional): number of different echoes that are used to calculate separate images.

These classifiers are used in the following notation (without spaces):

<Magnetisation modification> – <Dimensionality of data acquisition> – <PULSE SEQUENCE type> – <Echo number>

4.2 PULSE SEQUENCE type

4.2.1 General

PULSE SEQUENCES create MR images whose signal behaviour is primarily determined by either a SPIN-ECHO (SE) or a GRADIENT-ECHO (GR). The PULSE SEQUENCE type is defined at the centre

of the 2D or 3D K-SPACE: If the K-SPACE centre is acquired with a SPIN-ECHO, then the PULSE SEQUENCE is classified as a SPIN-ECHO PULSE SEQUENCE in this International Standard. If the centre of K-SPACE is not acquired with a SPIN-ECHO, then the PULSE SEQUENCE is classified as a GRADIENT-ECHO PULSE SEQUENCE in this International Standard.

A more precise classification the PULSE SEQUENCE type is achieved by providing information about the other K-SPACE lines. Therefore, the number of SPIN-ECHOES and GRADIENT-ECHOES per RF excitation are given.

Some image characteristics are sensitive to whether K-SPACE is totally acquired after a single RF excitation (single shot) or whether multiple RF excitations are required (multi-shot), so the number of RF excitations (shots) is also given.

4.2.2 Notation

For a SPIN-ECHO PULSE SEQUENCE the following notation is used:

<PULSE SEQUENCE type> : (SE_<Index1> GR_<Index2>)_<Index3>

For a GRADIENT-ECHO PULSE SEQUENCE the following notation is used:

<PULSE SEQUENCE type> : (GR_<Index1> SE_<Index2>)_<Index3>

Here, <Index1> and <Index2> denote the number of SPIN-ECHOES and GRADIENT-ECHOES, respectively. The values of <Index1> and <Index2> can be either integer numbers or formulas of the variables N and M using the signs “+”, “-”, “/”, and “x”. Here, N describes the total number of K-SPACE lines, and M is the number of K-SPACE lines per excitation. The part “GR_<Index2>” or “SE_<Index2>” is omitted, if <Index2> is zero.

<Index3> denotes the number of required RF excitations and is given either as an integer number or a formula of the variables N and M using the signs “+”, “-”, “/”, and “x”.

NOTE The sum of <Index1> and <Index2> is typically known as the echo train length, and <Index3> is often written as the number of shots.

4.3 Magnetisation modification

4.3.1 General

Optionally, the imaging characteristics of a PULSE SEQUENCE can be changed by adding gradients and RF pulses or by replacing parts of the PULSE SEQUENCE. These added pulses and gradients are used to prepare the magnetisation.

Different magnetisation modifications can be combined. The magnetisation modification is either acting primarily on the LONGITUDINAL MAGNETISATION or the TRANSVERSE MAGNETISATION, or is done during RF excitation.

For classification of the magnetisation modification, all parts of the PULSE SEQUENCES that have an additional influence on the image characteristics and that are not included in the classification of the PULSE SEQUENCE type are listed.

4.3.2 Notation

The magnetisation modification is written as a series of symbols as given in Table 1:

<Magnetisation modification> : <Symbol1> – <Symbol2> – ... – <SymbolN>

If no magnetisation modification is used, this classifier is omitted.

Table 1 – Magnetisation modification techniques

Symbol	Name	Physical principle
IR	Inversion recovery	Inversion of the LONGITUDINAL MAGNETISATION
SR	Saturation recovery	Saturation of the LONGITUDINAL MAGNETISATION
T2P	T2 preparation	Storage of a T2 contrast in the LONGITUDINAL MAGNETISATION using e.g. the pulses (90°)-(180°)-(-90°)
T2SP	T2* preparation	Storage of a T2* contrast in the LONGITUDINAL MAGNETISATION using e.g. the pulses (90°)-TE-(-90°)
SSAT	Spectral (chemical) saturation	Spectrally selective saturation of spins (e.g. fat, silicone, water)
SIR	Spectral (chemical) inversion	Spectrally selective inversion of spins (e.g. fat, silicone, water)
MTC	Magnetisation transfer contrast	Indirect saturation of the water signal using magnetisation transfer of spins bound to macromolecules
RSAT	(Regional) pre-saturation, tagging	Spatially selective saturation
RLAB	(Regional) labelling, spin labelling	Spatially selective excitation or inversion for tagging of moving spins
DE	Driven equilibrium	RF pulse at the end of the echo train for restoring the LONGITUDINAL MAGNETISATION
NS	Non-spatially selective excitation	RF pulse without slice selection gradient to excite all magnetisation in the volume of the transmitting RF coil
NRRF	Profiled (non-rectangular) RF pulse	RF excitation pulse with a dedicated, non-rectangular slice profile (e.g. to avoid saturation effects in TOF MR angiographies or for simultaneous excitation of two slices)
2DRF	2D selective RF pulse	RF pulse with a spatial selectivity in 2 dimensions
SSRF	Spatial-spectral RF pulse	Simultaneous spatially and spectrally selective excitation
DIFF	Diffusion weighting	Signal attenuation by additional gradients that induce a signal dephasing for diffusing spins
FLOWCn	Flow compensation ^a	Compensation of the n th gradient moment to suppress flow related signal changes, where n is a positive integer
FLOWSn	Flow sensitisation	Amplification of the n th gradient moment to enhance flow related signal changes, where n is a positive integer
T1R	Spin lock	T1 ρ contrast by additional RF pulses
SPOIL	TRANSVERSE MAGNETISATION spoiling	Spoiling by gradient and/or radio frequency
AREF	Refocusing of all gradients	Refocusing of all gradients within one TR interval (balanced steady state free precession)
PREF	Partial refocusing of the gradients	Refocusing of some gradients within one TR interval
NREF	Refocusing of the gradients in the next TR interval	Refocusing of the gradients so that an echo is acquired in the next TR interval
OFFSET	Offset echo	Time offset between the SPIN ECHO and acquisition of the centre of K-SPACE

^a If the classifier FLOWCn is not explicitly given, the PULSE SEQUENCE does not utilise flow compensation.

4.4 Dimensionality

4.4.1 General

With three-dimensional data acquisitions the spatial encoding of the signal from a volume is performed in all three dimensions, and data are sorted into a 3D K-SPACE. In two-dimensional data acquisitions encoding is performed only in the two in-plane dimensions of a selected slice, and a 2D K-SPACE is filled.

4.4.2 Notation

The dimensionality for a two-dimensional data acquisition is given by

<Dimensionality> : 2D

For a three-dimensional data acquisition the following notation is used:

<Dimensionality> : 3D

4.5 Echo number

4.5.1 General

Following an RF excitation, data can be acquired at different echo times, which leads to MR images with different contrasts. The number of reconstructed images with different contrasts is classified by the echo number.

4.5.2 Notation

The following notation is used for the echo number:

<Echo number> : E_ <Index>

<Index> is an integer number. If <Index> has a value of 1, the classifier is omitted.

Annex A (informative)

Examples of use of the PULSE SEQUENCE classification

Table A.1 shows examples of the application of this International Standard to MANUFACTURER-specific PULSE SEQUENCE names.

Table A.1 – MANUFACTURER-specific classification examples

GE	Hitachi	Philips	Siemens	Toshiba	This International Standard
EPI	GE EPI	Single shot FFE-EPI	FID-EPI	EPI	2D-(GR_N)_1
EPI	SE EPI	Single shot SE-EPI	SE-EPI	EPI	2D-(SE_1-GR_N-1)_1
FGRE-ET	Multi shot EPI	Multi shot EPI	Segmented EPI	Multi-shot EPI	2D-(GR_M)_N/M
SPIN-ECHO	SPIN-ECHO	SPIN-ECHO	SPIN-ECHO	SPIN-ECHO	2D-(SE_1)_N
SS-FSE	Single shot FSE	Single shot TSE	HASTE/RARE	FASE	2D-(SE_N)_1
fast SPIN-ECHO	fast SPIN-ECHO	turbo SPIN-ECHO	turbo SPIN-ECHO	FSE	2D-(SE_M)_N/M
dual echo fast SPIN ECHO		<td>	<td>		2D-(SE_M)_N/M-E_2
FLAIR/STIR	FLAIR/STIR	FLAIR/STIR	FLAIR/STIR	FLAIR/STIR	IR-2D-(SE_M)_N/M
GRASE	-	GRASE	TGSE	Hybrid EPI	2D-(SE_1-GR_M-1)_N/M
Spoiled GRASS	RSSG	T1-FFE	FLASH	FE	SPOIL-2D-(GR_1)_N
GRASS	SARGE	FFE	FISP	FE	PREF-2D-(GR_1)_N
FIESTA	BASG	Balanced FFE	True FISP	True SSFP	AREF-2D-(GR_1)_N
SSFP	TRSG	T2-FFE	PSIF	SSFP	NREF-2D-(GR_1)_N

Bibliography

DICOM 2008 standard: *Digital Imaging and Communications in Medicine (DICOM)* [viewed 2010-07-15]. Available from <ftp://medical.nema.org/medical/dicom/2008/>

RadLex Term Browser [viewed 2010-07-15]. Available from <http://www.radlex.org/viewer>

Index of defined terms used in this standard

NOTE The defined terms used in this International Standard may be looked up at <http://std.iec.ch/glossary>.

GRADIENT-ECHO (GR) 3.6
K-SPACE 3.4
LONGITUDINAL MAGNETISATION 3.3
MAGNETIC RESONANCE (MR) IEC 60601-2-33:2010, 201.3.217
MANUFACTURER IEC 60601-1:2005, 3.55
MR EQUIPMENT IEC 60601-2-33:2010, 201.3.218
PULSE SEQUENCE 3.1
SPIN-ECHO (SE) 3.5
TRANSVERSE MAGNETISATION 3.2

British Standards Institution (BSI)

BSI is the national body responsible for preparing British Standards and other standards-related publications, information and services.

BSI is incorporated by Royal Charter. British Standards and other standardization products are published by BSI Standards Limited.

About us

We bring together business, industry, government, consumers, innovators and others to shape their combined experience and expertise into standards-based solutions.

The knowledge embodied in our standards has been carefully assembled in a dependable format and refined through our open consultation process. Organizations of all sizes and across all sectors choose standards to help them achieve their goals.

Information on standards

We can provide you with the knowledge that your organization needs to succeed. Find out more about British Standards by visiting our website at bsigroup.com/standards or contacting our Customer Services team or Knowledge Centre.

Buying standards

You can buy and download PDF versions of BSI publications, including British and adopted European and international standards, through our website at bsigroup.com/shop, where hard copies can also be purchased.

If you need international and foreign standards from other Standards Development Organizations, hard copies can be ordered from our Customer Services team.

Subscriptions

Our range of subscription services are designed to make using standards easier for you. For further information on our subscription products go to bsigroup.com/subscriptions.

With **British Standards Online (BSOL)** you'll have instant access to over 55,000 British and adopted European and international standards from your desktop. It's available 24/7 and is refreshed daily so you'll always be up to date.

You can keep in touch with standards developments and receive substantial discounts on the purchase price of standards, both in single copy and subscription format, by becoming a **BSI Subscribing Member**.

PLUS is an updating service exclusive to BSI Subscribing Members. You will automatically receive the latest hard copy of your standards when they're revised or replaced.

To find out more about becoming a BSI Subscribing Member and the benefits of membership, please visit bsigroup.com/shop.

With a **Multi-User Network Licence (MUNL)** you are able to host standards publications on your intranet. Licences can cover as few or as many users as you wish. With updates supplied as soon as they're available, you can be sure your documentation is current. For further information, email bsmusales@bsigroup.com.

BSI Group Headquarters

389 Chiswick High Road London W4 4AL UK

Revisions

Our British Standards and other publications are updated by amendment or revision.

We continually improve the quality of our products and services to benefit your business. If you find an inaccuracy or ambiguity within a British Standard or other BSI publication please inform the Knowledge Centre.

Copyright

All the data, software and documentation set out in all British Standards and other BSI publications are the property of and copyrighted by BSI, or some person or entity that owns copyright in the information used (such as the international standardization bodies) and has formally licensed such information to BSI for commercial publication and use. Except as permitted under the Copyright, Designs and Patents Act 1988 no extract may be reproduced, stored in a retrieval system or transmitted in any form or by any means – electronic, photocopying, recording or otherwise – without prior written permission from BSI. Details and advice can be obtained from the Copyright & Licensing Department.

Useful Contacts:

Customer Services

Tel: +44 845 086 9001

Email (orders): orders@bsigroup.com

Email (enquiries): cservices@bsigroup.com

Subscriptions

Tel: +44 845 086 9001

Email: subscriptions@bsigroup.com

Knowledge Centre

Tel: +44 20 8996 7004

Email: knowledgecentre@bsigroup.com

Copyright & Licensing

Tel: +44 20 8996 7070

Email: copyright@bsigroup.com



...making excellence a habit.™