



BSI Standards Publication

Soil quality — Measurement of enzyme activity patterns in soil samples using fluorogenic substrates in micro-well plates

National foreword

This Draft for Development is the UK implementation of ISO/TS 22939:2010.

This publication is not to be regarded as a British Standard.

It is being issued in the Draft for Development series of publications and is of a provisional nature. It should be applied on this provisional basis, so that information and experience of its practical application can be obtained.

Comments arising from the use of this Draft for Development are requested so that UK experience can be reported to the international organization responsible for its conversion to an international standard. A review of this publication will be initiated not later than 3 years after its publication by the international organization so that a decision can be taken on its status. Notification of the start of the review period will be made in an announcement in the appropriate issue of *Update Standards*.

According to the replies received by the end of the review period, the responsible BSI Committee will decide whether to support the conversion into an international Standard, to extend the life of the Technical Specification or to withdraw it. Comments should be sent to the Secretary of the responsible BSI Technical Committee at British Standards House, 389 Chiswick High Road, London W4 4AL.

The UK participation in its preparation was entrusted to Technical Committee EH/4, Soil quality.

A list of organizations represented on this committee can be obtained on request to its secretary.

This publication does not purport to include all the necessary provisions of a contract. Users are responsible for its correct application.

© BSI 2010

ISBN 978 0 580 66288 1

ICS 13.080.30

Compliance with a British Standard cannot confer immunity from legal obligations.

This Draft for Development was published under the authority of the Standards Policy and Strategy Committee on 31 October 2010.

Amendments issued since publication

| Date | Text affected |
|------|---------------|
|------|---------------|

TECHNICAL SPECIFICATION

DD ISO/TS 22939:2010

ISO/TS 22939

First edition
2010-03-15

Soil quality — Measurement of enzyme activity patterns in soil samples using fluorogenic substrates in micro-well plates

*Qualité du sol — Mesure en microplaques de l'activité enzymatique
dans des échantillons de sol en utilisant des substrats fluorogènes*



Reference number
ISO/TS 22939:2010(E)

© ISO 2010

PDF disclaimer

This PDF file may contain embedded typefaces. In accordance with Adobe's licensing policy, this file may be printed or viewed but shall not be edited unless the typefaces which are embedded are licensed to and installed on the computer performing the editing. In downloading this file, parties accept therein the responsibility of not infringing Adobe's licensing policy. The ISO Central Secretariat accepts no liability in this area.

Adobe is a trademark of Adobe Systems Incorporated.

Details of the software products used to create this PDF file can be found in the General Info relative to the file; the PDF-creation parameters were optimized for printing. Every care has been taken to ensure that the file is suitable for use by ISO member bodies. In the unlikely event that a problem relating to it is found, please inform the Central Secretariat at the address given below.



COPYRIGHT PROTECTED DOCUMENT

© ISO 2010

All rights reserved. Unless otherwise specified, no part of this publication may be reproduced or utilized in any form or by any means, electronic or mechanical, including photocopying and microfilm, without permission in writing from either ISO at the address below or ISO's member body in the country of the requester.

ISO copyright office
Case postale 56 • CH-1211 Geneva 20
Tel. + 41 22 749 01 11
Fax + 41 22 749 09 47
E-mail copyright@iso.org
Web www.iso.org

Published in Switzerland

Contents

Page

| | |
|---|-----------|
| Foreword | iv |
| Introduction | v |
| 1 Scope | 1 |
| 2 Normative references | 1 |
| 3 Abbreviated terms | 1 |
| 4 Principle | 1 |
| 5 Reagents | 2 |
| 6 Apparatus and materials | 5 |
| 7 Procedure | 6 |
| 8 Calculation of results | 7 |
| 9 Expression of results | 8 |
| 10 Test report | 8 |
| Annex A (informative) Guidance on the use of freshly prepared substrates | 9 |
| Annex B (informative) Example of a graph for calculation | 11 |
| Bibliography | 13 |

Foreword

ISO (the International Organization for Standardization) is a worldwide federation of national standards bodies (ISO member bodies). The work of preparing International Standards is normally carried out through ISO technical committees. Each member body interested in a subject for which a technical committee has been established has the right to be represented on that committee. International organizations, governmental and non-governmental, in liaison with ISO, also take part in the work. ISO collaborates closely with the International Electrotechnical Commission (IEC) on all matters of electrotechnical standardization.

International Standards are drafted in accordance with the rules given in the ISO/IEC Directives, Part 2.

The main task of technical committees is to prepare International Standards. Draft International Standards adopted by the technical committees are circulated to the member bodies for voting. Publication as an International Standard requires approval by at least 75 % of the member bodies casting a vote.

In other circumstances, particularly when there is an urgent market requirement for such documents, a technical committee may decide to publish other types of document:

- an ISO Publicly Available Specification (ISO/PAS) represents an agreement between technical experts in an ISO working group and is accepted for publication if it is approved by more than 50 % of the members of the parent committee casting a vote;
- an ISO Technical Specification (ISO/TS) represents an agreement between the members of a technical committee and is accepted for publication if it is approved by 2/3 of the members of the committee casting a vote.

An ISO/PAS or ISO/TS is reviewed after three years in order to decide whether it will be confirmed for a further three years, revised to become an International Standard, or withdrawn. If the ISO/PAS or ISO/TS is confirmed, it is reviewed again after a further three years, at which time it must either be transformed into an International Standard or be withdrawn.

Attention is drawn to the possibility that some of the elements of this document may be the subject of patent rights. ISO shall not be held responsible for identifying any or all such patent rights.

ISO/TS 22939 was prepared by Technical Committee ISO/TC 190, *Soil quality*, Subcommittee SC 4, *Biological methods*.

Introduction

Micro-organisms are responsible for many key processes in the cycle of elements. Enzymes are responsible for the degradation of organic molecules and their mineralization. The main postulate is the microbial origin of soil enzymes, even if plant root exudates include enzymes. Extracellular enzymes in soil play key roles in the biodegradation of organic macromolecules. The simultaneous monitoring of several enzyme activities important in the biodegradation of organic compounds and mineralization of C, N, P and S in soil may reveal harmful effects caused by chemicals and other anthropogenic impacts. However, the measurements carried out under selected laboratory conditions using artificial substrates cannot be a substitute for the actual rate of enzymatic processes in soil *in situ*.

Soil quality — Measurement of enzyme activity patterns in soil samples using fluorogenic substrates in micro-well plates

1 Scope

This Technical Specification specifies a method for the measurement of several enzyme activities simultaneously in soil samples. Enzyme activities of soil vary seasonally and depend on the chemical, physical and biological characteristics of soil. Its application for the detection of harmful effects of toxic chemicals or other anthropogenic impacts depends on the simultaneous comparison of enzyme activities in a control soil similar to the test soil, or on exposure tests with chemicals or treatments.

2 Normative references

The following referenced documents are indispensable for the application of this document. For dated references, only the edition cited applies. For undated references, the latest edition of the referenced document (including any amendments) applies.

ISO 10381-6, *Soil quality — Sampling — Part 6: Guidance on the collection, handling and storage of soil under aerobic conditions for the assessment of microbiological processes, biomass and diversity in the laboratory*

ISO 10390, *Soil quality — Determination of pH*

ISO 10694, *Soil quality — Determination of organic carbon and total carbon after dry combustion (elementary analysis)*

3 Abbreviated terms

| | |
|------|--|
| E.C. | Enzyme code number defined by the Nomenclature Committee of the International Union of Biochemistry and Molecular Biology (NC-IUBMB) |
| SOM | Soil organic matter content |
| MUB | Modified universal buffer |

4 Principle

This Technical Specification describes a method for the simultaneous measurements of several enzymes in soil samples. It is based on the use of soil samples diluted in buffer containing fluorogenic substrates, which are incubated for 3 h at $(30 \pm 2)^\circ\text{C}$ in multi-well plates. After the incubation the enzyme activities are measured as fluorescence with a plate-reading fluorometer (References [1] and [2] in the Bibliography). The method described is based on dried standard and substrate plates enabling storage and limiting bias due to differences between reagent batches, and also enabling comparison between reagent batches. Annex A describes a method utilizing freshly prepared reagents, which has a clearly defined and exact incubation period. The advantage of the use of freshly prepared substrates is that an instrument for lyophilization is not required.

5 Reagents

5.1 Buffers

5.1.1 General

The selection of the buffer depends on the soil sample because the pH strongly affects enzyme activities. Sodium acetate buffer, 0,5 mol/l, at pH 5,5 has been used for acid soils with a high organic matter content. The use of the modified universal buffer (MUB) at the pH of the soil sample gives the flexibility necessary for coverage of a broad spectrum of different soils. Adequate stability of substrates at different buffers needs to be ensured. Good stability has been observed in 0,5 mol/l sodium acetate buffer at pH 5,5 (Reference [3] in the Bibliography).

5.1.2 Sodium acetate buffer, 0,5 mol/l, pH 5,5

- sodium acetate trihydrate 68,04 g;
- deionized water, ad 1 000 ml;
- acetic acid > 99,8 %.

Dissolve sodium acetate trihydrate in water (e.g. 800 ml) and adjust the pH to 5,5 with concentrated acetic acid (> 99,8 %; pro-analysis). Fill up to 1 000 ml. Sterilize in an autoclave at $(121 \pm 3) ^\circ\text{C}$ for 20 min. Store in a refrigerator for a maximum of two weeks.

5.1.3 Modified universal buffer (MUB) (Reference [4])

5.1.3.1 Stock solution

- tris(hydroxymethyl)aminomethane 12,1 g;
- maleic acid 11,6 g;
- citric acid 14,0 g;
- boric acid 6,3 g;
- sodium hydroxide (1 mol/l) 488 ml;
- deionized water, ad 1 000 ml.

Dissolve the ingredients and store the solution in a refrigerator.

5.1.3.2 Final buffer

- hydrochloric acid (0,1 mol/l);
- sodium hydroxide (0,1 mol/l).

Place 200 ml of the stock solution (5.1.3.1) in a 500 ml beaker containing a magnetic stirring bar, and place the beaker on a magnetic stirrer. Set the required pH with hydrochloric acid or with sodium hydroxide. Adjust the volume to 1 000 ml with deionized water. Sterilize in an autoclave at $(121 \pm 3) ^\circ\text{C}$ for 20 min.

5.2 Substrates and standards

5.2.1 Preparation of standard solutions

5.2.1.1 4-Methylumbelliferone (MUF) solution

- 4-methylumbelliferone (MUF) 0,022 0 g;
- dimethylsulfoxide (DMSO), ad 25 ml.

MUF in powder form can be stored at room temperature but protected from light. Weigh MUF carefully and dissolve it in DMSO in a brown volumetric flask, avoiding exposure to daylight. The solution cannot be stored.

5.2.1.2 7-Amino-4-methylcoumarin (AMC) solution

- 7-amino-4-methylcoumarin (AMC) 0,021 9 g;
- dimethylsulfoxide (DMSO), ad 25 ml.

AMC as powder can be stored in the refrigerator. Weigh AMC carefully and dissolve it in DMSO in a brown volumetric flask, avoiding exposure to daylight. The solution cannot be stored.

5.2.2 Preparation of substrate solutions

Commercially available fluorogenic substrates are delivered as powders that can be stored deep-frozen at (-20 ± 2) °C. On the day of use, weigh the amount required for a 1 000 $\mu\text{mol/l}$, 2 500 $\mu\text{mol/l}$ or 2 750 $\mu\text{mol/l}$ concentration in a volume of, for example, 50 ml, avoiding exposure to light. Weigh the powder into a brown volumetric flask and fill to the required volume with DMSO.

The volume should be big enough for reliable weighing and measurement of volumes. It also depends on the number of plates needed.

The commonly used dispensers are able to distribute simultaneously just one volume (e.g. 40 μl) to eight rows. To facilitate the use of these instruments enabling good volumetric precision, 2 500 $\mu\text{mol/l}$ solutions of the substrates should be prepared. However, for 4-MUF- β -D-glucopyranoside and for 4-MUF-phosphate substrates, a solution with the concentration of 2 750 $\mu\text{mol/l}$ is needed in order to produce the same final concentration of 500 $\mu\text{mol/l}$. These two solutions are further diluted simultaneously with the addition of the sample; 20 μl dimethylsulfoxide is added to the wells of these two substrates to facilitate dissolution. For chitinase activity measurement, a lower concentration is needed in order to avoid substrate inhibition, and the preparation of a solution with a concentration of 1 000 $\mu\text{mol/l}$ 4-MUF-*N*-acetyl- β -D-glucosaminide can be used to produce the final concentration of 200 $\mu\text{mol/l}$.

5.2.3 Preparation of multi-well plates

The substrate and standard solutions are added to multi-well plates as solutions and dried (e.g. freeze-dried) on the multi-well plates directly after dispensing. Dry plates can be stored at (-20 ± 2) °C for a year. Exposure to light shall be avoided during handling and storage of substrates and standards. A separate multi-well plate for substrates and standards has proved to be convenient.

5.2.4 Preparation of standard plates

Adequate replicate measurements, e.g. three to four replicates, are necessary due to the small sample volume. Standardization requires several concentrations of MUF or AMC, in replicate. Exposure to light shall be avoided during the dilution of standards. Calculate the required volume that depends on the number of samples and multi-well plates prepared. One example for the preparation of standards covering a wide range of enzyme activities is given below, but modifications can be made depending on the range of enzyme activities in the samples studied.

The stock solution of MUF with a concentration of 5 mmol/l is used to produce the dilutions containing 1 000 µmol/l, 500 µmol/l, 250 µmol/l, 125 µmol/l, 50 µmol/l, 25 µmol/l and 5 µmol/l MUF. Distribute the volumes needed (e.g. 40 ml) into a multi-well plate for concentrations of 0 µmol/l, 1,0 µmol/l, 5,0 µmol/l, 10 µmol/l, 25 µmol/l, 50 µmol/l, 100 µmol/l and 200 µmol/l, in replicate. This step is critical for the measurement uncertainty.

NOTE 1 This set of stock solutions enables the use of automatic dispensers, which yield a significantly better precision than manual pipetting.

The stock solution of AMC with a concentration of 5 mmol/l is used to produce the dilutions containing 250 µmol/l, 125 µmol/l, 50 µmol/l, 25 µmol/l, 5 µmol/l, 2,5 µmol/l and 0,5 µmol/l AMC. Distribute the volumes needed (e.g. 40 ml) into a multi-well plate for concentrations of 0 µmol/l, 0,1 µmol/l, 0,5 µmol/l, 1,0 µmol/l, 5,0 µmol/l, 10 µmol/l, 25 µmol/l and 50 µmol/l, in replicate. This step is critical for the measurement uncertainty.

NOTE 2 This set of stock solutions enables the use of automatic dispensers, which yield a significantly better precision than manual pipetting.

5.2.5 Preparation of substrate plates

Exposure to light shall be avoided during dilution of substrates. When using multi-well plates with dry substrates with the final substrate concentration of 500 µmol/l and the sample volume of 200 µl, a volume of 40 µl of the 2 500 µmol/l solution for the substrates is added, in replicate, to the wells. For 4-MUF-β-D-glucopyranoside and for 4-MUF-phosphate, a 2 750 µmol/l substrate solution is added, in replicate, to the wells. For chitinase activity measurement, 40 µl of the substrate 4-MUF-N-acetyl-β-D-glucosaminide is added as a 1 000 µmol/l solution to reach the final concentration of 200 µmol/l.

This is the concentration that has been used for several different soils with the assumption of an approximate saturation level. In validation tests for a broad spectrum of soils, appropriate substrate concentrations should be checked and/or an enzyme kinetic approach considered.

If a plate with 96 wells and an automatic dispenser are used for eight substrates, 12 replicates are conveniently available. When using four replicates, it is possible to analyse three different samples or dilution levels on one plate.

5.2.6 Drying of multi-well plates

Handling depends on the instrument e.g. lyophilizator used. Exposure to light shall be avoided during processing of the multi-well plates.

5.2.7 Fluorogenic substrates

Table 1 gives a list of fluorogenic substrates and standards that are available commercially¹⁾.

1) Glycosynth and Sigma are examples of producers of fluorogenic molecules. This information is given for the convenience of users of this document and does not constitute an endorsement by ISO of the producer named. Equivalent products may be used if they can be shown to lead to the same results.

**Table 1 — Fluorogenic artificial substrates available commercially
for the enzyme activity measurements**

| Enzyme | NC-IUBMB (see Clause 3) | Substrate | Element | Macromolecule degraded |
|---|----------------------------|---|------------|--|
| Arylsulfatase | E.C. 3.1.6.1 | 4-MUF-sulfate | Sulfur | Mineralization of organic sulfur |
| α -Glucosidase | E.C. 3.2.1.20 | 4-MUF- α -D-glucopyranoside | Carbon | Starch and glycogen |
| Cellobiosidase | E.C. 3.2.1.91 | 4-MUF- β -cellobiopyranoside | Carbon | Cellulose |
| β -Xylosidase | E.C. 3.2.1.37 | 4-MUF- β -D-xylopyranoside | Carbon | Xylane, xylobiose |
| β -Glucosidase | E.C. 3.2.1.21 | 4-MUF- β -D-glucopyranoside | Carbon | Cellulose |
| Phosphodiesterase (PDE) | E.C. 3.1.4.1 | bis-(4-MUF)-phosphate | Phosphorus | Hydrolysis of phosphate diesters |
| Chitinase | E.C. 3.2.1.30 | 4-MUF- <i>N</i> -acetyl- β -D-glucosaminide | Carbon | Breaking β -1-4-glycosidic bonds in <i>N</i> -acetyl-glucosaminide (chitin) and chitobiose |
| Phosphomonoesterase (PME) | E.C. 3.1.3.2 | 4-MUF-phosphate | Phosphorus | Hydrolysis of phosphate monoesters |
| Leucine-aminopeptidase | E.C. 3.4.11.1 | L-leucine-AMC | Nitrogen | Hydrolysis of oligopeptides → aminoacids |
| Alanine-aminopeptidase | E.C. 3.4.11.12 | L-alanine-AMC | Nitrogen | Hydrolysis of oligopeptides → aminoacids |
| MUF = 4-methylumbelliferone AMC = 7-amino-4-methylcoumarin | | | | |

6 Apparatus and materials

6.1 Equipment for the homogenization of the soil samples.

6.1.1 Sieves, with grid size, e.g. 4 mm.

NOTE Other grid sizes can be used depending on the soil texture.

6.1.2 Mechanical homogenizer²⁾.

For some soil types, an ultrasonic disaggregator has been used.

6.2 Usual laboratory glassware.

6.3 Multi-well plates³⁾, with covers.

2) Bamix rod homogenizer and OmniMixer are examples of suitable homogenizers. This information is given for the convenience of users of this document and does not constitute an endorsement by ISO of these products. Equivalent products may be used if they can be shown to lead to the same results.

3) There are differences in background fluorescence between brands and MUF and AMC analyses may require different brands.

6.4 Automatic dispenser for reagents (optional)⁴⁾.

NOTE Compared with manual pipetting, an automatic dispenser decreases significantly the uncertainty of volumes dispensed.

6.5 Instrument for drying of the multi-well plates (e.g. by lyophilization).

6.6 Incubators, set at (30 ± 2) °C and allowing the use of a plate shaker at 450 min^{-1} to 700 min^{-1} .

6.7 Plate-reading fluorometer, with the excitation wavelength set at 355 nm and the emission wavelength set at 460 nm, with the excitation lamp energy set appropriately⁵⁾.

7 Procedure

7.1 Sampling

Take and handle soil samples as specified in ISO 10381-6. A composite sample consisting of 20 sub-samples taken in the field and homogenized by sieving through the selected sieves has been observed to yield reasonably low uncertainty of measurement for soil samples.

NOTE This method has been applied to boreal forest soil samples with a high soil organic matter (SOM) content and to agricultural soils with a high clay or silt content. For coarse samples with a relatively big particle size, a different sample handling might be necessary.

Soil pH is an important characteristic affecting enzyme activities and shall be measured in accordance with ISO 10390. This information can be used in the evaluation of enzyme activity results, even if the measurements are not carried out at *in situ* pH. The soil organic-matter content correlates strongly with soil enzyme activities. It shall be measured in accordance with ISO 10694.

Storage in a refrigerator after sieving for two days is not suitable but samples can be stored deep-frozen at (-20 ± 2) °C for at least four months.

7.2 Sample preparation

7.2.1 Homogenization

The sieved soil is homogenized, and a test sample of 4 g is added to 120 ml of the selected buffer at (22 ± 2) °C. Treatment with the homogenizer (6.1.2) for 3 min in an ice bath and treatment with the homogenizer using $9\,600 \text{ min}^{-1}$ for 3 min in an ice bath have been used. The homogenized sample is adjusted into a final volume of 200 ml (dilution 1:50).

For some soil types, an ultrasonic disaggregator has been used: disperse the soil sample in the buffer (dilution dependent on the activity of the soil) and homogenize the soil suspension using an ultrasonic disaggregator for 120 s at an output energy of $50 \text{ J}\cdot\text{s}^{-1}$ (Reference [5] in the Bibliography).

4) Wallac 1298-003 Delfia is an example of a suitable instrument. This information is given for the convenience of users of this document and does not constitute an endorsement by ISO of this instrument. Equivalent instruments may be used if they can be shown to lead to the same results.

5) Wallac Victor is an example of a suitable instrument allowing excitation lamp energy to be set optimally. This information is given for the convenience of users of this document and does not constitute an endorsement by ISO of this instrument. Equivalent instruments may be used if they can be shown to lead to the same results.

7.2.2 Preparation of dilutions

The optimal dilution level depends on the soil sample and enzyme. Dilutions 1:100 and 1:1 000 are usually adequate. It is advisable to use the same dilution level for each enzyme because different dilutions do not give exactly comparable results. This may be due to the differences in exposure of enzymes to substrates in microsites.

When using plates with dried substrates and standards.

- Prepare the dilution of 1:100 by adding 20 ml of tempered buffer (5.1.2 or 5.1.3) at (22 ± 2) °C to 20 ml of homogenized sample (7.2.1).
- Prepare the dilution of 1:1 000 by adding 36 ml of tempered buffer (5.1.2 or 5.1.3) at (22 ± 2) °C to 4 ml of the dilution 1:100 prepared above.

7.2.3 Sample distribution

In four replicates, add 200 μ l of the diluted soil sample to each well of the multi-well plate containing freeze-dried substrates, to yield substrate concentrations of 500 μ mol/l. Add 20 μ l of DMSO to the wells containing substrates for β -glucosidase and phosphomonoesterase. Place the cover on each plate.

Add 200 μ l in three or, preferably, four replicates of the homogenized, diluted sample to the wells with each standard. Place the cover on each plate.

A separate standard curve is necessary for each soil sample and dilution for MUF and AMC separately.

The blank values are obtained by measuring the substrate plates immediately after adding the sample.

NOTE This kind of blank does not reveal chemical instability of fluorogenic compounds. This is the reason why the stability of the compounds in sterile buffer needs to be tested separately.

7.3 Incubation

Incubate the multi-well plates for 3 h at (30 ± 2) °C while shaking continuously, e.g. at about 700 min^{-1} .

NOTE The incubation temperature affects reaction rates and the optima depends on the enzymes. Depending on the study, a different temperature, as described in this Technical Specification, can be used, e.g. the *in situ* temperature.

7.4 Fluorescence measurements

Measure the fluorescence directly after adding the sample to the multi-wells and after a 3 h incubation, using excitation at 355 nm and emission at 460 nm with adjusted excitation lamp energy.

NOTE Fluorogenic compounds dried in the multi-wells can yield elevated fluorescence directly after addition of the sample to the micro-wells because of undissolved compounds on the well surfaces. In soil samples, this is usually not a significant source of measurement uncertainty. However, it is also possible to measure fluorescence repeatedly to yield information on the dissolution, reaction rates and kinetic parameters (Reference [6] in the Bibliography).

The stability of the substrates shall be tested for each substrate batch (powder from the manufacturer) by measuring the change in fluorescence in the sterile buffer during incubation (Reference [1] in the Bibliography).

8 Calculation of results

The standard curve is plotted for MUF or AMC molar concentration (μ mol/l) versus fluorescence. The MUF or AMC concentration of blank (c_b) and sample (c_s) are read from the standard curve.

An example of a graph is given in Annex B.

The result is calculated by subtracting the average of four measurement replicates of blanks from the sample and multiplying the difference by the dilution factor, using either soil volume, soil fresh mass, dry mass or soil organic matter (SOM).

9 Expression of results

The results are expressed as micromoles per litre ($\mu\text{mol/l}$) of MUF or AMC released per soil volume or per gram (g) of soil fresh mass and/or dry mass and/or soil organic matter during 3 h, depending on the study. The measurements necessary for each expression, volumetric or gravimetric measurements or the determination of loss on ignition for soil organic matter shall be carried out.

Soil characteristics vary widely due to geography, climate and land use. The interpretation of results cannot, currently, be based on set limit values for each enzyme activity. The experimental design shall facilitate comparisons with a control soil or between samples from relevant sites.

10 Test report

The test report should include the following information:

- a) a reference to this Technical Specification;
- b) adequate identification of the sample;
- c) details of storage temperature and duration;
- d) type of soil and soil physical and chemical characterization;
- e) pH value of the soil sample;
- f) buffer and incubation conditions applied;
- g) test results;
- h) any details not specified in this Technical Specification or which are optional, as well as any incident which may have influenced the results.

Annex A (informative)

Guidance on the use of freshly prepared substrates

A.1 Introduction

The fluorogenic substrates can also be applied as freshly prepared substrates according to Reference [7] in the Bibliography. It is preferable to analyse fresh soil samples directly after sampling and homogenization. However, because the reagents are not stable in solution and variation between batches may be significant, it is sometimes advisable to store all the soil samples to be compared as weighted deep-frozen aliquots, and to carry out the measurements of enzyme activities on all the soil samples to be compared simultaneously, rather than to carry out repeated measurements using separate batches of reagents. The weighing of substrate and standard compounds is a potential source of measurement uncertainty and an increase in the mass to be weighed decreases the measurement uncertainty.

A.2 Reagents

A.2.1 Buffer

The following buffers have been applied for the measurement of enzyme activities using freshly prepared substrates.

The selection of the buffer is critical. The buffering capacity shall be adequate to enable the set pH to be constant throughout the measurement. Either a constant pH or a pH optimized for the soil in question, or optimized separately for each enzyme in question, can be selected. However, the feasibility of enzyme activity pattern analysis supports the use of the same buffer for all the enzyme activities measured simultaneously. The same buffers as those described in 5.1.2. and 5.1.3 can also be used for freshly prepared substrates and the soil sample can be diluted directly in the buffer as in 7.2.1. MES (2-[*N*-morpholino]ethanesulfonic acid) buffer at pH 6,1 for the majority and Trizma buffer⁶⁾ at pH 7,8 for the aminopeptidase activities are commonly used with freshly prepared substrates.

Use MES buffer for phosphatase and enzymes involved in C cycling (substrates containing MUF as the fluorescent compound, see Table 1). To obtain 0,1 mol/l MES buffer (pH 6,1), dissolve 22,1 g in 1 l of water.

Use Trizma buffer for peptidases (substrates containing AMC as the fluorescent compound, see Table 1). To obtain 0,05 mol/l Trizma buffer (pH 7,8), dissolve 0,985 g of Trizma base and 2,66 g of Trizma HCl in 0,5 l of water.

Sterilize in an autoclave at (121 ± 3) °C for 20 min.

6) Trizma buffer is an example of a suitable product available commercially. This information is given for the convenience of users of this document and does not constitute an endorsement by ISO of this product.

A.2.2 Substrates

Stock solution (10 mmol/l): dissolve substrates in 300 µl of dimethylsulfoxide (DMSO) and adjust the volume to 10 ml with autoclaved water.

Working solutions (1 mmol/l): dilute stock solution with autoclaved MES buffer for MUF substrates and autoclaved Trizma buffer for AMC substrates. The total volume depends on the number of samples. For one sample series, all the substrate solutions should be produced at once.

NOTE When the sample volume is equal to the substrate volume in each well, the final substrate concentration is 500 µmol/l.

A.2.3 Standards

Dissolve standards in dimethylsulfoxide to concentrations of 5 mmol/l; subsequently dilute with MES buffer (for substrates containing MUF as the fluorescent compound) or Trizma buffer (for substrates containing AMC as the fluorescent compound) to a final concentration of 10 µmol/l. The total volume depends on the number of samples. For one sample series, all the standard solutions should be produced at once.

A.3 Procedure

A.3.1 Substrate plates

The substrate solutions can be added to multi-well plates as solutions, following the addition of the sample dilution to the wells.

Disperse 1 g of soil with 100 ml of sterile deionized water (for optimal dilution, see 7.2.2) using an ultrasonic disaggregator (50 J·s⁻¹ for 120 s). Mix 50 µl of the soil suspension with 50 µl of the appropriate autoclaved buffer (MES buffer for MUF substrates, Trizma for AMC substrates) and 100 µl of substrate solution in microplates. Three or four replicates of each sample is preferable.

NOTE Manual pipetting of 50 µl volumes of substrates, and especially pipetting of soil samples in suspension in 50 µl aliquots, causes very high measurement uncertainty.

A.3.2 Standard plates

The standard solutions can be added to multi-well plates as solutions following the addition of the sample dilution to the wells.

Mix standards with 50 µl of soil suspension (separate standard curve for each sample, see 7.2.3) and the appropriate amount of buffer to obtain final concentrations of 0 µmol/l, 0,5 µmol/l, 1 µmol/l, 2,5 µmol/l, 4 µmol/l and 6 µmol/l in a final volume of 200 µl. To account for quenching, mix 100 µl of buffer with 100 µl of substrate for each substrate.

A.3.3 Incubation

For incubation and fluorescence measurements, see 7.3 and 7.4.

When freshly prepared plates are used with the fluorogenic compound already in solution, even a very short incubation period can yield measurable enzyme activities and the change in fluorescence shall be measured in precise constant time intervals.

Annex B (informative)


Example of a graph for calculation

The data for standardization for MUF and AMC are given in Table B.1. The graphs drawn and curves fitted using a second class polynomial for MUF are represented in Figure B.1, and for AMC in Figure B.2. Other curve-fitting approaches may be needed depending on the sample and fluorogen.

NOTE It is sometimes possible to pool data for the standardization curve from different samples if soil characteristics are similar.


Table B.1 — Fluorescence measurement results for MUF and AMC standards

| MUF | | AMC | |
|------------------------------------|---------|------------------------------------|---------|
| Concentration $\mu\text{mol/g}$ | Counts | Concentration $\mu\text{mol/g}$ | Counts |
| 0 | 9 652 | 0 | 2 807 |
| 0 | 9 837 | 0 | 2 771 |
| 0 | 9 629 | 0 | 2 918 |
| 1 | 12 309 | 0,1 | 3 185 |
| 1 | 13 054 | 0,1 | 3 188 |
| 1 | 12 532 | 0,1 | 3 223 |
| 5 | 29 430 | 0,5 | 4 224 |
| 5 | 27 881 | 0,5 | 4 496 |
| 5 | 29 940 | 0,5 | 4 375 |
| 10 | 50 703 | 1 | 7 091 |
| 10 | 53 686 | 1 | 6 671 |
| 10 | 54 361 | 1 | 6 892 |
| 25 | 112 601 | 5 | 44 442 |
| 25 | 128 289 | 5 | 45 015 |
| 25 | 128 303 | 5 | 42 234 |
| 50 | 240 218 | 10 | 94 225 |
| 50 | 241 312 | 10 | 92 720 |
| 50 | 239 976 | 10 | 95 688 |
| 100 | 342 290 | 25 | 223 674 |
| 100 | 421 423 | 25 | 218 396 |
| 100 | 416 506 | 25 | 237 820 |
| 200 | 676 064 | 50 | 306 714 |
| 200 | 692 148 | 50 | 329 967 |
| 200 | 665 945 | 50 | 335 326 |



Key
X concentration of MUF, in $\mu\text{mol/g}$
Y fluorescence signal, in counts

Figure B.1 — Standardization curve for MUF



Key
X concentration of AMC, in $\mu\text{mol/g}$
Y fluorescence signal, in counts

Figure B.2 — Standardization curve for AMC

Bibliography

- [1] VEPSÄLÄINEN, M., KUKKONEN, S., VESTBERG, M., SIRVIÖ, H., NIEMI, R.M., 2001. Application of soil enzyme activity test kit in a field experiment. *Soil Biology and Biochemistry* **33**, pp. 1665-1672
- [2] VEPSÄLÄINEN, M., ERKOMAA, K., KUKKONEN, S., VESTBERG, M., NIEMI, R.M., 2004. The impact of crop plant cultivation and peat amendment on soil microbial activity and structure. *Plant and Soil* **264**(1-2), pp. 273-286
- [3] NIEMI, R.M., VEPSÄLÄINEN, M., 2005. Stability of the fluorogenic enzyme substrates and pH optima of enzyme activities in different Finnish soils. *Journal of Microbiological Methods* **60** (2), pp. 195-205
- [4] TABATABAI, M.A., 1994. *Soil enzymes*. In: *Methods of Soil Analysis. Part 2. Microbiological and Biochemical Properties*. R.W. Weaver, S. Angle, P. Bottomley, D. Bezdicek, S. Smith, A. Tabatabai and A. Wollum, editors, pp. 775-833. Soil Science Society of America, Madison
- [5] STEMMER, M., GERZABEK, M.H., KANDELER, E., 1998. Organic matter and enzyme activity in particle-size fractions of soils obtained after low-energy sonication. *Soil Biology and Biochemistry* **30**, pp. 9-17
- [6] MARX, M.-C., KANDELER, E., WOOD, M., WERMBTER, N., JARVIS, S.C., 2005. Exploring the enzymatic landscape: distribution and kinetics of hydrolytic enzymes in soil particle-size fractions. *Soil Biology and Biochemistry* **37**, pp. 35-48
- [7] MARX, M.-C., WOOD, M., JARVIS, S.C., 2001. A microplate fluorimetric assay for the study of enzyme diversity in soils. *Soil Biology and Biochemistry* **33**, pp. 1633-1640

ICS 13.080.30

Price based on 13 pages

British Standards Institution (BSI)

BSI is the national body responsible for preparing British Standards and other standards-related publications, information and services.

BSI is incorporated by Royal Charter. British Standards and other standardization products are published by BSI Standards Limited.

About us

We bring together business, industry, government, consumers, innovators and others to shape their combined experience and expertise into standards-based solutions.

The knowledge embodied in our standards has been carefully assembled in a dependable format and refined through our open consultation process. Organizations of all sizes and across all sectors choose standards to help them achieve their goals.

Information on standards

We can provide you with the knowledge that your organization needs to succeed. Find out more about British Standards by visiting our website at bsigroup.com/standards or contacting our Customer Services team or Knowledge Centre.

Buying standards

You can buy and download PDF versions of BSI publications, including British and adopted European and international standards, through our website at bsigroup.com/shop, where hard copies can also be purchased.

If you need international and foreign standards from other Standards Development Organizations, hard copies can be ordered from our Customer Services team.

Subscriptions

Our range of subscription services are designed to make using standards easier for you. For further information on our subscription products go to bsigroup.com/subscriptions.

With **British Standards Online (BSOL)** you'll have instant access to over 55,000 British and adopted European and international standards from your desktop. It's available 24/7 and is refreshed daily so you'll always be up to date.

You can keep in touch with standards developments and receive substantial discounts on the purchase price of standards, both in single copy and subscription format, by becoming a **BSI Subscribing Member**.

PLUS is an updating service exclusive to BSI Subscribing Members. You will automatically receive the latest hard copy of your standards when they're revised or replaced.

To find out more about becoming a BSI Subscribing Member and the benefits of membership, please visit bsigroup.com/shop.

With a **Multi-User Network Licence (MUNL)** you are able to host standards publications on your intranet. Licences can cover as few or as many users as you wish. With updates supplied as soon as they're available, you can be sure your documentation is current. For further information, email bsmusales@bsigroup.com.

BSI Group Headquarters

389 Chiswick High Road London W4 4AL UK



Revisions

Our British Standards and other publications are updated by amendment or revision.

We continually improve the quality of our products and services to benefit your business. If you find an inaccuracy or ambiguity within a British Standard or other BSI publication please inform the Knowledge Centre.

Copyright

All the data, software and documentation set out in all British Standards and other BSI publications are the property of and copyrighted by BSI, or some person or entity that owns copyright in the information used (such as the international standardization bodies) and has formally licensed such information to BSI for commercial publication and use. Except as permitted under the Copyright, Designs and Patents Act 1988 no extract may be reproduced, stored in a retrieval system or transmitted in any form or by any means – electronic, photocopying, recording or otherwise – without prior written permission from BSI. Details and advice can be obtained from the Copyright & Licensing Department.

Useful Contacts:

Customer Services

Tel: +44 845 086 9001

Email (orders): orders@bsigroup.com

Email (enquiries): cservices@bsigroup.com

Subscriptions

Tel: +44 845 086 9001

Email: subscriptions@bsigroup.com

Knowledge Centre

Tel: +44 20 8996 7004

Email: knowledgecentre@bsigroup.com

Copyright & Licensing

Tel: +44 20 8996 7070

Email: copyright@bsigroup.com

...making excellence a habit.™