



INTERNATIONAL STANDARD ISO 4074:2002 TECHNICAL CORRIGENDUM 2

Published 2008-04-15

Corrected version 2002-12-01

INTERNATIONAL ORGANIZATION FOR STANDARDIZATION • МЕЖДУНАРОДНАЯ ОРГАНИЗАЦИЯ ПО СТАНДАРТИЗАЦИИ • ORGANISATION INTERNATIONALE DE NORMALISATION

Natural latex rubber condoms — Requirements and test methods

TECHNICAL CORRIGENDUM 2

Préservatifs masculins en latex de caoutchouc naturel — Exigences et méthodes d'essai

RECTIFICATIF TECHNIQUE 2

Technical Corrigendum 2 to ISO 4074:2002 was prepared by Technical Committee ISO/TC 157, *Mechanical contraceptives*.

Introduction

This corrigendum shall be read together with ISO 4074:2002 and Technical Corrigendum 1 (2003).

Page 4, 6.2

Delete the whole paragraph.

Page 10, Table A.1

In column 1, row 3 delete the term “oven-treated”

Page 11, Table B.1

In column 1, row 3 delete the term “oven-treated”

Page 22, H.1

The clause now reads:

Oven treatment is used to condition condoms for shelf-life determination. This annex describes the method for oven treatment.