



**INTERNATIONAL STANDARD ISO 6009:1992**  
**TECHNICAL CORRIGENDUM 1**

Published 2008-03-01

INTERNATIONAL ORGANIZATION FOR STANDARDIZATION • МЕЖДУНАРОДНАЯ ОРГАНИЗАЦИЯ ПО СТАНДАРТИЗАЦИИ • ORGANISATION INTERNATIONALE DE NORMALISATION

## **Hypodermic needles for single use — Colour coding for identification**

### **TECHNICAL CORRIGENDUM 1**

*Aiguilles hypodermiques non réutilisables — Code de couleurs pour l'identification*

*RECTIFICATIF TECHNIQUE 1*

Technical Corrigendum 1 to ISO 6009:1992 was prepared by Technical Committee ISO/TC 84, *Devices for administration of medicinal products and intravascular catheters*.

---

*Page 2*

Replace Annex A with the following.

**Annex A**  
(informative)

**Reference colour samples**

The reference colour samples referred to in this International Standard are available by contacting:

Melab  
Medizintechnik und Labor GmbH  
Mollenbachstrasse 19  
71229 Leonberg  
Germany  
info@melab.de

This information is given for the convenience of the users of this International Standard and does not constitute an endorsement by ISO of the reference colour samples.